





L.E.D. X Ray Viewer Box (Economical Version)



Bio-X specializes in the manufacture of X-ray film viewers (X-ray illuminator).


Keeping the needs of the doctors in mind along with the continuously changing demands of hospitals and Nursing homes, we have developed an economical slim and effective L.E.D. Viewing system to view M.R.I, C.T. and X Ray images. The economical version of Bio-X slim LED X ray viewer is available in Single Screen and Double Screen.

We offer a combination of benefits that is unmatched by our competition: reliability, distinctive product line designed to match your requirements, technical support for layouts, excellent service and stunning looks.

This extremely slim model has the power adaptor outside the view box.

FEATURES

- Slim, Durable and Effective
- Energy Efficient: Savings range from 82% to 93%
- No UV emission
- Unparallel latest LED technology 4-5 times longer than regular CCFL/EEFL
- Power efficient. L.E.D Light with long life- upto 100,000 hours
- 8 steps digital dimmer to get the perfect brightness
- External power adapter



Parameter		Single Screen	Double Screen
Dimensions (mm)	Frame (L x H x T)	430(W)X490(H)X16(D)mm	763(W)X490(H)X16(D)mm
	Viewing Area (W x H)	360 x 440	736 x 440
Power Consumption		15W	30W
Light Source		LED	LED
Power Supply Input		AC 100 - 240V 50/ 60 Hz	AC 100 - 240V 50/ 60 Hz
Brightness (LUX)		> 6,500	> 6,500
Net Weight (KGS)		3.6	7