


A petit product gives you a perfect performance in your professional profile, Enjoy It




Bldg.46-1, BDA International Business Park, Beijing100176, CHINA

Tel.: +86-10-67856941 / 67856340 Ext-840 Fax: +86-10-67856343

<http://www.vhmedical.com> email: [sales@vhmedical.com](mailto:sales@vhmedical.com)

ISO13485



# iCardioView

**Innovative, Intuitive and Simple**

**A unique gadget meets for all your professional demands in cardiac world**

**Simultaneous 12-Lead ECG**

**Automatic Measurements and Interpretations**


**Filters: Low-Pass (different setting)**

**Baseline Wander**

**Motion Artifact**

**Report: Blue Tooth Printer or PDF upload**

**Download:**



## ECG Recorder Specification

- Sampling Rate:  
A/D: 24K SPS/Ch  
Recording: 1K SPS/Ch  
Patented Multi-Channel Synchronous A/D
- Quantization Precision:  
A/D: 24 Bits  
Recording: 16 Bits
- Resolution: 0.4 $\mu$ V
- Common Mode Rejection: >110dB
- Input Impedance: >20M $\Omega$
- Frequency Response: 0.05-250Hz( $\pm$ 3dB)
- Time Constant:  $\geq$ 3.2Sec.
- Maximum Electrode Potential:  $\pm$ 300mV DC
- Dynamic Range:  $\pm$ 15mV
- Defibrillation Protect: Build-in
- Communication: WiFi
- Power: 2 $\times$ AA Batteries

